


Sample Preparation Tools (WRB 1401, assisted)

Mirror Suggested for the Clones Observation

SP Bel-Art Petri Dish Plate Reader:

- This reader allows hands-free viewing of Plates or Petri dishes through a large, adjustable mirror with 3× magnification.
- Applications include counting of cell cultures for colony or plaque growth. The reader assists in evaluation of monolayer cell growth in Petri dishes, viewing of color changes, and locating cells in microtiter plates (VWR Cat#12574-064)



Bullet (Cluster) Tubes Suggested for Small Volume Samples on the FCRC Analyzers

Corning® 96-Well Polypropylene Tubes:

- Composed of 96 PP tubes in a standard microplate format
- 1.2 ml tubes available individually or in strips of 8 tubes (VWR Cat#29442-604, bulk; #29442-612, sterile in racks)
- PE tube caps available in 8-cap strips (#29442-614)

Aspirator with Multi-Channel Adapter

Corning™ 8-Channel Adapter for Vacuum Aspirator:

- Perfect for safe removal/disposal of supernatants from a wide range of laboratory vessels (Corning Cat#4930)
- For use with Corning Vacuum Aspirator. This 8-Channel Adapter can be used to remove media and solutions from 8 wells simultaneously on 96 well plate (Corning Cat#4931)

